# Orthocem

**■ DENTSCARE LTDA** - Av. Edgar Nelson Meister, 474 - Distrito Industrial - 89219-501 - Joinville - SC - Authorization of Functioning MS PSX44XYOXX28 - CNPJ/Tax ID: 05.106.945/0001-06 - BRAZILIAN INDUSTRY - Registration at ANVISA Nº 80172310043 Technical Authority: Friedrich Georg Mittelstädt - CRQ: 13100147-SC - Brand: FGM®

Rev: 02

Irritante  
Irritant  
Irritant

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## **EN Instructions Manual**

---

### **Cement/Adhesive for Orthodontic Brackets Fixation Professional Use Only**

Read all this information carefully before using the product. Keep it for later reference until total consumption of the product and or until there is no more interaction between the product and your last patient.

#### **Product Description**

Orthocem is an orthodontic single component cement/adhesive that promotes the union between metal and/or ceramics brackets to the enamel surface. The product is lightcuring, proportioning adequate work time to adjust bracket's position. Another characteristic is that the product contains fluoride, preventing dental demineralization process, and it also has an ideal viscosity to facilitate the brackets positioning, preventing it's displacement before the fotovotación.

The product's application form also has a clinical advantage: primer and bond are united in the syringe, allowing the bracket's clinical installation to be resumed to enamels acid etch (etch and rinse) and adhesive/cement application on the bracket, followed by light activation.

#### **Presentation Form**

Package containing 1 syringe of 4g, 1 phosphoric acid (CondAc 37 - FGM) - Registro na ANVISA nº 80172310025 and instruction's manual.

#### **Basic Composition**

Contains methacrylic monomers (BisGMA, TEGDMA) and methacrylic phosphatized monomers, stabilizers, cantoroquinone, co-initiators and nanometric load of silicon

dioxide.

#### **Product Indication**

Fixation of metal, ceramic or polycarbonate brackets to enamel surface.

#### **Precautions and Contraindications**

- Professional Use Only;
- Avoid contact to eyes and skin. Wear gloves to manipulate the product;
- In case of contact with the eyes, wash immediately with water in abundance. If there is persistent irritation, search for medical support;
- Wear protecting goggles during the light activation of the product;
- In case of contact with skin, rinse with water and soap. If necessary, search medical support;
- Avoid contact of phosphoric acid to mucosa, skin and eyes. It can cause burns. Wear gloves for its manipulation;
- Do not expose the material to high temperatures or intense light.

#### **Side Effects**

Do not store the product near eugenol-based materials because those can inhibit the adequate cure of adhesive/cement.

#### **Instructions**

Before initiating the treatment, read carefully the Precautions, Contraindications and possible Side Effects.

1. Isolation of operatory field: the operatory field must be kept free of saliva contamination, by means of an efficient isolation. Additionally, it's recommended using ArcFlex (FGM) lip retractor for better access to the operatory field.
2. Teeth preparation: proceed to teeth prophylaxis with oil free pastes. Rinse with water in abundance and dry the teeth with air free of oil and humidity.
3. Acid etch: do the acid etch on enamel with 37% phosphoric acid (CondAc 37 - FGM) following manufacturer's instructions.
4. Rinse with water in abundance and air dry.
5. Bracket bonding: apply little amount of Orthocem straight from the syringe onto the base of the bracket and put it on position on previously etched enamel.
6. Remove the excess of adhesive without displacing it from its correct position.
7. For bonding brackets made of metal, ceramic or polycarbonate, lightcure the adhesive for 40 seconds in each bracket varying between the margins, with the lightcuring unit as near as possible to the adhesive.
- NOTE: Considering a photo-polymerization with a minimum power of 450 mW/cm<sup>2</sup> and blue light with 400-500nm.

8. The arches can be placed immediately after bracket bonding.

9. When concluding the orthodontic treatment, the brackets should be removed with the aid of specific instruments. The residual cement adhered to dental surface can be removed with low rotation specific burs or sandpaper disks of fine grain, according to the case. The removal of the mentioned residues should be done carefully to avoid damages to sound enamel.

10. Then, the teeth should be polished with felt disks (Diamond Flex - FGM) and polishing paste of extra-fine grain (Diamond Excel - FGM).

#### **Conservation and Storage**

Do not expose the product to high temperatures or to light. Do not store the material next to other containing eugenol because it can inhibit the adequate polymerization process of the adhesive. Store the product at a temperature between 5 a 30°C. Do not use the product after its expiry date.

#### **Warnings**

Do not use the product after its expiry date. For disposal follow your country's legislation. Do not reuse the empty package. Keep out of children's reach. Avoid the contact with skin and eyes. It is important to use procedure gloves when manipulating the material.

*This product was manufactured for dental use only and it should be manipulated according to the instructions. The manufacturer is not responsible for damage caused by other types of use or by*

*incorrect manipulation. Furthermore, the user is obligated to verify, before the utilization and under his/her responsibility, if this product is compatible with the desired utilization, especially when such use is not mentioned in the instructions. Data descriptions do not constitute any type of warranty and, therefore, are not legally binding.*